

- MT1X-
- MTnnnn-1X-
- MTnnnnX-1X- (but not MTN02X)

**ProM2<sup>®</sup> human MHC Class II Monomer:** In order to identify antigen-specific CD4<sup>+</sup> T lymphocytes, fluorochrome-labeled Class II tetramers are required. ProM2<sup>®</sup> human MHC Class II Monomer reagents can be made into Class II tetramers when combined with Streptavidin fluorochrome conjugates. Streptavidin has four biotin-binding sites, enabling biotinylated ProM2<sup>®</sup> human MHC Class II Monomer reagents to form Class II tetramers. CD4<sup>+</sup> T cells stained with Class II tetramers can be analyzed by flow cytometry and the frequency of antigen-specific T cells determined.

**For Research Use Only. Not for use in therapeutic or diagnostic procedures.**

**Test size:** ProM2<sup>®</sup> Monomers are provided in 35 µg and 100 µg sizes.

**Concentration/Formulation:** ProM2<sup>®</sup> human MHC Class II Monomer is supplied at 0.4 mg/ml in 20 mM Tris, pH 8, 50 mM NaCl, stabilized with 0.5% BSA and 0.025% sodium azide.

**Storage Condition:** -80°C. **Avoid freeze-thaw cycles.**

**Shelf Life:** Use liquid nitrogen to flash-freeze upon receipt of material, or tetramerize immediately. The monomer is stable for 12 months if stored as instructed above.

**Hazards:** This reagent is formulated in 0.025% sodium azide. Under acid conditions the toxic compound hydrazoic acid may be released. Compounds containing sodium azide should be flushed with running water while being discarded.

#### Quality Control Assay Results

**Appearance:** Colorless solution

**Protein Characterization:** Passed

**Released by:**  
(Date as per product label above)

#### Class II Tetramer Production Protocol

**Additional materials required:** Streptavidin-R-PE or Streptavidin-APC for ProM2<sup>®</sup> Monomer, PBS containing 0.025% sodium azide.

1. Spin Streptavidin-R-PE or Streptavidin-APC in a chilled microcentrifuge at 14,000 ×g for 3 minutes. This will remove protein aggregates that contribute to non-specific staining. Maintain reagents on ice, shielded from light, until required. Do not aspirate any part of the pelleted aggregates when taking test volumes for conjugation.
2. *To conjugate 35 µg ProM2<sup>®</sup> human MHC Class II Monomer with Streptavidin-R-PE:*  
Add 13 µl of 0.8 mg/ml Streptavidin-R-PE to 35 µg ProM2<sup>®</sup> human MHC Class II Monomer, mix gently and incubate at 4°C for 15 minutes. Repeat the addition of Streptavidin-R-PE four times with a 15 minute gap between each addition. Make up to a final volume of 400 µl with PBS/0.025% sodium azide.
3. *To conjugate 100 µg ProM2<sup>®</sup> human MHC Class II Monomer with Streptavidin-APC:*  
Add 23 µl of 0.09 mg/ml Streptavidin-APC to 35 µg ProM2<sup>®</sup> human MHC Class II Monomer, mix gently and incubate at 4°C for 15 minutes. Repeat the addition of Streptavidin-APC four times with a 15 minute gap between each addition. Make up to a final volume of 400 µl with PBS/0.025% sodium azide.
3. *To conjugate 100 µg ProM2<sup>®</sup> human MHC Class II Monomer with Streptavidin-R-PE:*  
Add 36 µl of 0.8 mg/ml Streptavidin-R-PE to 100 µg ProM2<sup>®</sup> human MHC Class II Monomer, mix gently and incubate at 4°C for 15 minutes. Repeat the addition of Streptavidin-R-PE four times with a 15 minute gap between each addition. Make up to a final volume of 1.15 ml with PBS/0.025% sodium azide.

To conjugate 100 µg ProM2<sup>®</sup> human MHC Class II Monomer with Streptavidin-APC:

Add 65 µl of 0.09 mg/ml Streptavidin-APC to 100 µg ProM2<sup>®</sup> human MHC Class II Monomer, mix gently and incubate at 4°C for 15 minutes. Repeat the addition of Streptavidin-APC four times with a 15 minute gap between each addition. Make up to a final volume of 1.15 ml with PBS/0.025% sodium azide.

Store Class II tetramers at 4°C, protected from light. **Do not freeze.**

We recommend that Class II tetramers are used at 5 µl / test.

### **Cellular Staining Protocol:**

**Additional materials required:** Wash buffer (0.1% sodium azide, 0.1% BSA in PBS), Fix solution (1% fetal calf serum, 2.5% formaldehyde in PBS), anti-CD4 antibody, anti-CD19 antibody (optional<sup>†</sup>).

1. Centrifuge Class II tetramer in a chilled microcentrifuge at 14,000 ×g for 5 minutes. This will remove protein aggregates that contribute to non-specific staining. Maintain reagents on ice, shielded from light, until required. Do not aspirate any part of the pelleted aggregates when taking tests for staining.
2. Allocate 1-2 × 10<sup>6</sup> lymphoid cells (PBMC) per staining condition.
3. Wash the cells with wash buffer and resuspend them in the residual volume (~50 µl).
4. Add 5 µl of Class II tetramer to the cells and mix by pipetting.
5. Incubate at 37°C for 2 hours in the dark.
6. Wash the cells in wash buffer.
7. Add an optimally titrated amount of anti-CD4 antibody per staining condition. At this stage, it is also recommended to add anti-CD19 antibody in order to gate out B cells when performing analysis.
8. Incubate on ice for 20-30 minutes in the dark.
9. Wash the cells twice in wash buffer and store them in fix solution in the dark until analysis.

The tetramer-positive cells are most conveniently viewed by gating first on live lymphoid cells and then analyzing on a two-color plot showing CD4 on the x-axis and tetramer on the y-axis.

<sup>†</sup>Tetramers can bind non-specifically to B cells. It is therefore strongly recommended to include anti-CD19 antibody when staining in order to gate on CD19<sup>-</sup> cells before plotting tetramer versus CD4.

### **Protocol Optimization:**

The following guidelines will help you optimize your protocol for the best possible results:

**Setting the live lymphocyte gate** It is important to ensure that the forward-scatter (FSC) and side-scatter (SSC) gates are set correctly on the cell population of interest. This is to ensure that dead cells, cell aggregates and cell debris are excluded from the fluorescence data.


**Titrating the Class II tetramer** Although a single test quantity of Class II tetramer should normally be sufficient to stain 1-2 × 10<sup>6</sup> cells, it is important that you first titrate the Class II tetramer. Carry out a range of doubling dilutions from 1 test per 1 × 10<sup>6</sup> cells down to 1/16 test per 1 × 10<sup>6</sup> cells.

**CD4 antibody** Investigate the effect of titrating the anti-CD4 antibody.

**Temperature** The temperature at which cells are stained can affect signal considerably. Varying time and temperature of incubation is necessary to determine optimal signal to noise ratio depending upon the MHC/peptide combination and T cell receptor. We recommend incubation at 37°C for 2 hours in the first instance.

**Positive control** Class II tetramers should be tested against a specific T cell line (or clone). Be sure to use T cells that have not been recently stimulated as this has been shown to cause down-regulation of T cell receptors.

**Negative Control** To control for non-specific staining it is also useful to stain T cells with CLIP-Class II tetramer (refer to MTN02 for further information). Alternatively, staining T cells from unexposed individuals may be used when detecting T cell responses to a specific antigen.



1 × 10<sup>6</sup> cells were incubated with 1 test size R-PE-labeled Class II tetramer at 37°C for 2 hours. Non-specific staining was eliminated from the plot by gating on CD19<sup>-</sup> cells before plotting CD4 vs Class II tetramer.